

Thermophilic helicase-dependent DNA amplification using the IsoAmp™ SE experimental kit for rapid detection of *Streptococcus equi* subspecies *equi* in clinical samples

Sergey Artiushin,¹ Yanhong Tong, John Timoney, Bertrand Lemieux, Anne Schlegel, Huimin Kong

Abstract. A simple and portable assay for detection of *Streptococcus equi* subspecies *equi* has been developed based on amplification of *S. equi*-specific sequence using a thermophilic helicase-dependent reaction followed by visual detection of the amplicon in a disposable lateral flow cassette. An experimental kit (IsoAmp™ SE) was evaluated. Analytical sensitivity was 50 copies of *S. equi* genomic DNA per reaction. The IsoAmp SE assay had 100% specificity when applied to nasal swabs and washes. The assay was more sensitive than culture but less sensitive than nested polymerase chain reaction (PCR). The test requires neither expensive equipment nor extensive training of personnel, provides a practical alternative to culture or PCR assays for detection of *S. equi* in clinical samples, and expedites identification of atypical colonies of *S. equi* and *Streptococcus zooepidemicus* in the laboratory.

Key words: Diagnosis; horses; strangles; *Streptococcus equi*; thermophilic helicase-dependent amplification.

Introduction

Strangles caused by *Streptococcus equi* subspecies *equi* is the most prevalent contagious bacterial disease of horses in the world. Pyrexia, followed by nasal discharge and abscessation of the lymph nodes of head and neck, are characteristic symptoms.^{20,22} In some horses, abscesses are formed in body organs distant from the head and neck.¹⁶ Susceptible horses usually acquire infection through direct horse-to-horse contact and indirectly via water sources or feed contaminated by nasal discharges. Shedding of *S. equi* usually begins after a latent period of 4–16 days and continues for 3–6 weeks after the acute phase. Some animals harbor *S. equi* in the guttural pouch for months after clinical recovery and may become long-term carriers that intermittently shed the organism.^{8,14,15} There are over 9.5 million horses in the United States, some of which are frequently moved long distances between race tracks and other competitive events and for breeding, sales, and training. Thus, there is potential for transmission of infectious agents including *S. equi*. Identification of horses persistently or intermittently shedding the infectious agent is crucially important to control strangles.^{18,28}

Prevention and control of strangles is highly dependent on knowledge of the infection status of individual horses. Culture of nasal swabs and washes, or pus aspirated from abscesses, which is considered the “gold standard” in detection of *S. equi*, may take several days, and requires significant expertise since

specimens often contain other beta-hemolytic streptococci. Moreover, identification of atypical *S. equi* requires additional biochemical tests.²⁰ Polymerase chain reaction (PCR) assay based on *seM*, the gene for the antiphagocytic M protein of *S. equi*,²⁴ is more rapid, sensitive, and specific than culture,^{11,14,23} but is technically complex, requires costly instrumentation, and reagents are not widely available. Also, shipment of samples to qualified laboratories delays decisions concerning movement and transportation of horses and increases the cost of testing.

Development of a rapid, specific, and affordable assay for a wider range of users may substantially expedite identification of *S. equi* in specimens from horses and improve control of the disease. Such an assay could be based on thermophilic helicase-dependent DNA amplification (tHDA). As for PCR, amplification of target sequence in tHDA is achieved by numerous cycles of DNA denaturation, primer annealing, and primers extension. However, tHDA uses heat-stable helicase to unwind complementary DNA strands instead of heat, allowing isothermal DNA amplification without a thermocycler.²⁷ Further simplification of the assay could be

From the Gluck Equine Research Center, University of Kentucky, Lexington, KY (Artiushin, Timoney), and the BioHelix Corporation, Beverly, MA (Tong, Lemieux, Schlegel, Kong).

¹Corresponding Author: Sergey Artiushin, Gluck Equine Research Center, University of Kentucky, Lexington, KY 40546. scart11@email.uky.edu

Table 1. Oligonucleotides used in the current study for polymerase chain reaction (PCR) and thermophilic helicase-dependent DNA amplification (tHDA) for rapid detection of *Streptococcus equi* in clinical samples.

Primer/Probe	Target gene	Oligonucleotide sequence (5'→3')	GenBank no.
SzpF (PCR forward)	<i>szPSe</i>	CAGCAGCAGAAGCAGCAATC	U73163
SzpR (PCR reverse)	<i>szPSe</i>	GTGTAGAAGTCGCTAGCACC	
Spr7 (PCR forward)	<i>seM</i>	TGCATAAAGAAGTTCCTGTC	U73162
Spr6 (PCR reverse)	<i>seM</i>	GATTCGGTAAGAGCTTGACG	
Spr8 (nested PCR forward)	<i>seM</i>	CGAACTCTGAGGTTAGTCGT	
Pr1 (nested PCR reverse)	<i>seM</i>	CATACCTACTCCATCAGCA	
15F (tHDA forward)	<i>se18.9</i>	AGCAACTGTAAGTACAGGAACAGTAAGT	DQ068464
5R (tHDA reverse)	<i>se18.9</i>	Bio-TTCTTAGGCGCAAATATTTCCATGAGTACC	
tHDA probe	<i>se18.9</i>	TCAAGATATAATTATACCGGTTGG-FITC	

attained by using a disposable lateral flow cassette for detection of the amplicon.^{12,27} Such a platform has been recently developed and commercialized^a in a kit format and utilized for identification of toxigenic *Clostridium difficile* in stool samples,⁷ *Helicobacter pylori*,¹⁷ *Vibrio cholera*, *Bacillus anthracis*,²⁶ and *Human immunodeficiency virus-1*.²¹ The assay takes less than 2 hr and can be performed in a minimally equipped laboratory or office setting.

In the current study, the development and evaluation of a tHDA assay for detection of *S. equi* in specimens from horses is described. The assay is based on asymmetrical amplification of the 5' end of *Se18.9* specific to *S. equi*²⁵ followed by hybridization of the amplicon to the labeled probe and detection of the hybridization product in a disposable cassette.

Material and methods

Bacterial culture, DNA isolation, and sample preparation

Twenty isolates of *S. equi* and 35 isolates of *Streptococcus zooepidemicus* isolated over the last decade (2000–2010) in different geographical locations were grown at 37°C in Todd–Hewitt broth (THB) supplemented with 0.2% yeast extract or THB agar. Streptococcal DNA were isolated from 5–10 ml of stationary-phase cultures. Following centrifugation at 10,000 × *g* for 3 min, bacterial pellets were resuspended in 200 µl of phosphate buffered saline (PBS) supplemented with 200 units of mutanolysin^b and incubated for 1 hr at 37°C. The DNA was isolated using a commercial kit^c from bacterial lysates following the manufacturer's instructions. Concentrations of DNA were determined using a spectrophotometer.^d

Development of *S. equi*-specific thermophilic helicase-dependent DNA amplification

Sequence encoding the N-terminal region of mature *Se18.9*²⁵ of *S. equi* with least homology in *S. zooepidemicus* was used to design primers and probe using the Primer3 program.¹⁹ Program parameters were set to select primers with melting

temperature, GC, and size of amplicon as 65–71°C, 42–47%, and 80–120 base pairs (bp), respectively. The reverse primer, 5R (Table 1), was labeled with biotin so that the resulting amplicon would bind to streptavidin-conjugated color particles for visualization (Fig. 1). The tHDA probe (Table 1) was labeled with fluorescein isothiocyanate (FITC) as a tag to anchor the colored amplicon-probe hybrid to FITC-specific antibody on the T-line of the detection strip in the disposable lateral flow cassette^a (Fig. 1). The primers and probes used in the current study were proprietary.^c An IsoAmpTM SE experimental kit for detection of *S. equi* was prepared by complementation of a commercial tHDA kit^a with primers 15F, 5R, and tHDA probe (Table 1). An excess of 5R primer was added to generate single-stranded DNA to bind to the hybridization probe in the final reaction product.

Analytical sensitivity of the IsoAmp SE kit

Ten, 100 copies of genomic DNA of 3 isolates of *S. equi* were used in triplicate to determine the analytical sensitivity of the kit. Twenty µl of streptococcal genomic DNA were mixed in a 0.2-ml tube with 16.5 µl of buffer mix,^a 3.5 µl of enzyme mix,^a and 10 µl of water. Following 90-min incubation at 65°C, the reaction tubes containing products of reaction were transferred into disposable cassettes for detection of amplicon by lateral flow capture assay (Fig. 1).

Analytical specificity of the IsoAmp SE kit

The analytical specificity of the assay was tested using 50 copies of genomic DNA from 20 *S. equi* isolates and 10³ and 10⁴ copies of genomic DNA from 35 isolates of *S. zooepidemicus* per reaction. The tHDA assay and visualization of results were performed as described above.

Evaluation of the IsoAmp SE kit for detection of *S. equi* in clinical specimens

Samples for tHDA and PCR were prepared from 38 clinical specimens as follows. One ml of nasal wash or 0.5 ml of PBS-containing material released from nasal swabs were centrifuged




Figure 1. Schematic diagram of the thermophilic helicase-dependent amplification (tHDA) assay using the IsoAmp™ SE kit^a and disposable lateral flow cassettes.^a **A**, the tHDA reaction is performed with an excess of a biotin-labeled primer 5R to produce a single-stranded amplicon that can hybridize to a probe containing a 3'-fluorescein isothiocyanate (FITC) reporter. Streptavidin-coated latex beads visualize the hybrid. A lateral flow strip with 2 capture zones is used to detect this product as well as validate the performance of the lateral flow device. Anti-FITC antibodies in the test line (T) capture the DNA complex, and biotin-conjugated to the control line (C) captures beads that have not bound to single strand amplicon. **B**, disposable cassettes showed positive and negative results.

at $10,000 \times g$ for 3 min. DNA was extracted from the pellets as described above using a commercial kit.^c Twenty μl of DNA was used for the tHDA reaction as described above.

Culture and nested PCR were used as reference controls during evaluation of the helicase-dependent DNA amplification experimental kits. For culture, clinical specimens were plated on THB agar supplemented with horse blood and incubated in 5% (v/v) carbon dioxide in air at 37°C for 48 hr and inspected at 24 and 48 hr for presence of mucoid, hemolytic colonies typical of *S. equi*. Additional fermentation tests were performed for non-mucoid colonies. *Streptococcus equi* was identified by its inability to ferment ribose, sorbitol, trehalose, and lactose.

Nested PCR was performed as follows. Aliquots (2 μl) of samples were added to 30 μl of $1 \times$ Taq DNA polymerase buffer^f containing 1 unit of Taq DNA polymerase, 2 mM of magnesium chloride, 0.2 mM deoxyribonucleotide triphosphate, and 0.2 μM of primers Spr7 and Spr6 (Table 1). The PCR was performed on a thermocycler^g using a program involving denaturation at 94°C for 2 min followed by 30 cycles of denaturation at 92°C for 1 min, primer annealing at 57°C for 30 sec, and primer extension at 72°C for 1 min. Following completion of the first round, 2 μl were transferred to a fresh reaction mix with nested primers Spr8 and Pr1 (Table 1), and a second round of PCR was performed using the same conditions. The PCR products were analyzed by ethidium bromide agarose electrophoresis. Amplicons produced by nested PCR of samples negative in tHDA and culture were purified and sequenced to confirm the presence of *seM* sequence of *S. equi*. Sequencing was performed by a commercial facility.^c

Results of tHDA and culture were scored in comparison with nested PCR as a reference standard. Sensitivity, specificity, and diagnostic accuracy were calculated as ratios of

(true positive/true positive + false negative), (true negative/true negative + false positive), and (true negative/true negative + false negative), respectively. True positive, false positive, true negative, and false negative were specimens that were PCR+/Test+, PCR-/Test+, PCR-/Test-, and PCR+/Test-, respectively.

Identification of atypical colonies of *S. equi* and *S. zooepidemicus*

Portions of colonies of non-mucoid isolates of *S. equi* or mucoid isolates of *S. zooepidemicus* were transferred using toothpicks into 0.5-ml tubes containing 15 μl of water. Three randomly selected colonies of each isolate of *S. equi* and *S. zooepidemicus* were used. Following incubation at 100°C for 10 min and centrifugation at $10,000 \times g$ for 2 min, 10- and 5- μl aliquots of supernatants were tested by tHDA and PCR, respectively. Primers SzpF and SzpR universal for *S. equi* and *S. zooepidemicus*, and Spr6 and Spr7 specific for *S. equi*, were used for PCR (Table 1). The primers used in the present study are proprietary.^h The tHDA and PCR assays were performed as described above.

Results

Development of *S. equi*-specific thermophilic helicase-dependent DNA amplification

Primers 15F and 5R (Table 1) were selected, following comparison of 20 primers, as the most specific and providing the highest amplicon yield in tHDA reaction. No homology to amplicon sequence has been detected in genomes of other microorganisms following a search of available data banks. An experimental kit, IsoAmp SE, was manufactured using

Table 2. *Streptococcus equi* subsp. *equi* and *Streptococcus zooepidemicus* isolate DNA reactivity using the IsoAmp™ SE kit^a used for rapid detection of *S. equi* in clinical samples.

Isolate	No. of isolates	Copies of genomic DNA per reaction		
		50	10 ³	10 ⁴
<i>S. equi</i>	20	Positive	Positive	Positive
<i>S. zooepidemicus</i>	32	Negative	Negative	Negative
	3	Negative	Negative	Positive

these primers and probe. Figure 1 illustrates the workflow of the assay including tHDA and lateral flow detection of amplicons. Asymmetric amplification of the target sequence in the presence of an excess of primer 5R results in enrichment of the reaction with single-stranded, biotin-labeled DNA. Such DNA is effectively hybridized at incubation temperature to the probe forming a biotin/FITC-labeled complex. The hybrid moves along the detection strip of the cassette, binds with red latex beads coated with streptavidin, and is captured by anti-FITC antibodies in the test line. Streptavidin-coated beads that have not been bound by amplicon are captured on the biotin control line. When only the control line is detected, the assay is scored as a true negative (Fig. 1).

Analytical sensitivity and specificity of the IsoAmp SE kit

Ten, 100 copies of genomic DNA of 3 isolates of *S. equi* were used in triplicate to determine the analytical sensitivity of the kit. Results were strongly positive in all experiments with quantities of template equivalent to 50 and 100 copies of *S. equi* genomic DNA per reaction. Reducing the quantity to 25 copies led to inconsistent results that are probably due to minor variation in concentration of target or capture antibody on the detection strip. No responses were detected for 10 copies of genomic DNA per reaction.

The analytical specificity of the assay was tested using DNA from 20 isolates of *S. equi* and 35 isolates of *S. zooepidemicus*. The DNA of *S. equi* isolates gave positive responses when 50 copies of genomic DNA were used in the assay (Table 2). No reactivity in tHDA was detected for 10³ copies of genomic DNA for each of 35 isolates of *S. zooepidemicus*. However, DNA of 3 isolates of *S. zooepidemicus* gave positive responses at 10⁴ copies per reaction.

Evaluation of thermophilic helicase-dependent DNA amplification assay for detection of *S. equi* in clinical samples

Polymerase chain reaction, tHDA, and culture were compared for detection of *S. equi* in equine samples. *Streptococcus equi* DNA was detected by tHDA in 20 out of 38 samples that

Table 3. Comparison of thermophilic helicase-dependent amplification (tHDA) and culture for detection of *Streptococcus equi* in clinical samples from horses.*

Test	True positive	False negative	False positive	True negative	Percent sensitivity	Percent specificity
tHDA	20	3	0	15	87	100
Culture	17	6	0	15	74	100

*True positive = PCR+/Test+; false negative = PCR+/Test-; false positive = PCR-/Test+; true negative = PCR-/Test-.

were also positive by nested PCR. No tHDA-positive samples were identified among PCR-negative samples showing that tHDA is as specific as PCR for detection of *S. equi*. Nevertheless, detection of *S. equi* DNA by nested PCR in 3 specimens negative by culture and tHDA revealed a higher sensitivity of the nested PCR. Sequence of amplicons produced by nested PCR of these specimens showed identity to the corresponding part of the *seM* gene of *S. equi*. Culture revealed *S. equi* in 17 out of 20 tHDA-positive specimens, whereas *S. equi* was not detected in tHDA-negative specimens. Thermophilic helicase-dependent DNA amplification was therefore as specific as culture but more sensitive in identification of *S. equi* in clinical samples (Table 3). Other Gram-positive flora including *S. zooepidemicus* contaminated all except 3 of the *S. equi*-positive and -negative specimens, but this did not interfere with assay specificity.

Comparison of tHDA and culture using the nested PCR as a standard revealed that the tHDA assay had a sensitivity of 87%, a specificity of 100%, and a diagnostic accuracy of 83%. Culture was as specific as PCR and tHDA for detection of *S. equi* in specimens from horses but had inferior sensitivity (74%) and a diagnostic accuracy of 71% (Table 3).

Identification of morphologically atypical *S. equi* and *S. zooepidemicus*

Nucleic acids extracted from a portion of colonies of 3 non-mucoid *S. equi* and 5 mucoid *S. zooepidemicus* were analyzed by tHDA and PCR. The presence of streptococcal DNA from colonies was confirmed by positive reactions in PCR with SzpF and SzpR primers specific for both *S. equi* and *S. zooepidemicus* (Table 4). *Streptococcus equi* isolates alone were positive in tHDA as well as in PCR when primers Spr6 and Spr7 specific for *seM* were used.

Discussion

In the current study, a tHDA-based assay (IsoAmp SE) for detection of *S. equi* in specimens from horses was developed and evaluated. The test required neither expensive equipment nor extensive training of personnel and may be used for point-of-care testing. The *Se18.9* gene encoding an antiphagocytic

Table 4. IsoAmp™ SE kit^a and polymerase chain reaction (PCR) identification of *Streptococcus equi* and *Streptococcus zooepidemicus* in clinical samples with atypical colony morphology.

Isolates	Morphology	Reactivity		
		IsoAmp SE	<i>szPSe</i> PCR	<i>seM</i> PCR
<i>S. equi</i>				
Sample 1	Non-mucoid	Positive	Positive	Positive
Sample 2	Non-mucoid	Positive	Positive	Positive
Sample 3	Non-mucoid	Positive	Positive	Positive
<i>S. zooepidemicus</i>				
Sample 1	Mucoid	Negative	Positive	Negative
Sample 2	Mucoid	Negative	Positive	Negative
Sample 3	Mucoid	Negative	Positive	Negative
Sample 4	Mucoid	Negative	Positive	Negative
Sample 5	Mucoid	Negative	Positive	Negative

virulence factor of *S. equi* that binds the H control factor of complement is secreted abundantly by *S. equi*, but not by the closely related *S. zooepidemicus*, which is commonly found in the nasal cavity of horses,²⁵ and was used as a target for the tHDA assay. Comparison of sequences of *Se18.9* from 8 isolates of *S. equi* and 10 isolates of *S. zooepidemicus* showed that its 5' terminal half is highly conserved among *S. equi* isolates of different geographical origin, whereas the ortholog in *S. zooepidemicus* is divergent and shares only limited homology with *Se18.9* (data not shown). Primers and a probe selected from this sequence revealed high specificity to *S. equi* giving positive results in tHDA with all DNA from 20 randomly selected isolates of this microorganism. The DNA from 35 isolates of *S. zooepidemicus* collected during the last decade from different geographical regions showed negative results in tHDA when 10³ copies of genomic DNA per reaction were used. However, additional experiments using 10⁴ copies of genomic DNA per reaction resulted in positive responses for 3 out of 35 isolates of *S. zooepidemicus* creating the potential for false-positive results in cases of heavy contamination by some *S. zooepidemicus* isolates (Table 2).

The analytical sensitivity of the tHDA assay was 50 copies of genomic DNA of *S. equi* per reaction and lies between that for conventional PCR and nested PCR, as previously described for *S. equi*.^{3,14} The sensitivity reported for this platform (combination of tHDA and disposable cassette) and that previously developed⁹ for identification of *Staphylococcus aureus* were the same. However, this sensitivity level is not the limit for this platform. For instance, the analytical sensitivity of the assay for detection of *C. difficile* was 20 copies of genomic DNA per reaction.⁷

Comparison of the tHDA test with assays currently used for identification of *S. equi* reveals that the specificity of tHDA is similar to that of culture and PCR when applied to nasal swabs and washes. Sensitivity and diagnostic accuracy of

tHDA were 13% higher than those of culture. Although evaluation of the tHDA applied to culture-negative clinical specimens from horses with different degrees of contamination did not reveal false-negative results (Table 3), a conclusion about the probability of false-positive results will require future trials on additional clinical samples.

A second potential application of the tHDA assay is rapid identification of atypical forms of group C streptococci periodically isolated from equine specimens. Culture is widely used for definitive diagnosis of strangles and should be routinely used to confirm diagnosis based on PCR. According to this technique, samples from horses are inoculated onto horse blood agar, incubated 24–48 hr, and inspected for hemolytic mucoid colonies characteristic of *S. equi*. The other Lancefield group C streptococci, *S. zooepidemicus* and *Streptococcus dysgalactiae* subsp. *equisimilis*, which frequently present in clinical samples, also form hemolytic colonies, but are usually non-mucoid. Ability to ferment lactose, trehalose, and sorbitol is an additional descriptor.^{13,20,22} Nevertheless, the correct identification of *S. equi* and other streptococci could be hampered in the case of isolates with atypical colony morphology and fermentation pattern.¹⁰ Several molecular methods including sequencing of the intergenic region of the 16S-23S RNA genes,⁶ multiplex PCR,² real-time PCR,⁴ multilocus sequence typing,²⁹ pulsed field gel electrophoresis,⁵ and random amplification of polymorphic DNA PCR^{1,5} have been evaluated for differentiation of *S. equi* and other streptococci. Although any one of these methods could be used for resolving ambiguous cases, the absence of essential reagents, equipment, and expertise in a diagnostic laboratory exclude practical application of these techniques. The current study demonstrated that simple sample preparation followed by tHDA could differentiate atypical *S. equi* and *S. zooepidemicus* in less than 2 hr.

In summary, a specific and sensitive assay for detection of *S. equi* in equine specimens based on tHDA reaction followed by visualization of the test result in a disposable cassette is described. The ability to perform this assay in a simple laboratory or even point-of-care by personnel with basic training could substantially reduce the time for making decisions concerning movement and/or isolation of horses. The assay could also be used in clinical laboratories to expedite identification of morphologically and biochemically atypical isolates of *S. equi* and *S. zooepidemicus*.

Sources and manufacturers

- IsoAmp®, BioHelix® Express Strip (BESt™) cassette type I, IsoAmp® II Universal tHDA Kit; BioHelix Corp., Beverly, MA.
- Sigma-Aldrich, St. Louis, MO.
- ZR Genomic DNA II Kit™, Zymo Research Corp., Irvine, CA.
- NanoDrop, Thermo Scientific, Wilmington, DE.
- Operon Biotechnologies Inc., Huntsville, AL.
- Invitrogen Corp., Carlsbad, CA.
- Sprint, Thermo Hybaid, Franklin, MA.
- Integrated DNA Technologies Inc., Coralville, IA.

Declaration of conflicting interests

The author(s) declared no potential conflicts of interest with respect to the research, authorship, and/or publication of this article.

Funding

The U.S. Department of Agriculture (USDA SBIR 2008-00036) and the Kentucky Science and Engineering Foundation (KSEF-148-502-08-228) supported the current research.

References

- Al-Ghamdi GM, Kapur V, Ames TR, et al.: 2000, Use of repetitive sequence-based polymerase chain reaction for molecular epidemiologic analysis of *Streptococcus equi* subspecies *equi*. *Am J Vet Res* 61:699–705.
- Alber J, El-Sayed A, Lammler C, et al.: 2004, Multiplex polymerase chain reaction for identification and differentiation of *Streptococcus equi* subsp. *zooepidemicus* and *Streptococcus equi* subsp. *equi*. *J Vet Med B* 51:455–458.
- Artiushin S, Timoney JF: 1997, PCR for detection of *Streptococcus equi*. *Adv Exp Med Biol* 418:359–361.
- Båverud V, Johansson SK, Aspan A: 2007, Real-time PCR for detection and differentiation of *Streptococcus equi* subsp. *equi* and *Streptococcus equi* subsp. *zooepidemicus*. *Vet Microbiol* 124:219–229.
- Bert F, Branger C, Lambert-Zechovsky N: 1997, Pulsed-field gel electrophoresis is more discriminating than multilocus enzyme electrophoresis and random amplified polymorphic DNA analysis for typing pyogenic streptococci. *Curr Microbiol* 34:226–229.
- Chanter N, Collin N, Holmes N, et al.: 1997, Characterization of the Lancefield group C streptococcus 16S-23S RNA gene intergenic spacer and its potential for identification and sub-specific typing. *Epidemiol Infect* 118:125–135.
- Chow WHA, McCloskey C, Tong Y, et al.: 2008, Application of isothermal helicase-dependent amplification with a disposable detection device in a simple sensitive stool test for toxigenic *Clostridium difficile*. *J Mol Diagn* 10:452–458.
- Fintl C, Dixon P, Brazil T, et al.: 2000, Endoscopic and bacteriological findings in a chronic outbreak of strangles. *Vet Rec* 147:480–484.
- Goldmeyer J, Li H, McCormac M, et al.: 2008, Identification of *Staphylococcus aureus* and determination of methicillin resistance directly from positive blood cultures by isothermal amplification and a disposable detection device. *J Clin Microbiol* 46:1534–1536.
- Grant S, Efstratiou A, Chanter N: 1993, Laboratory diagnosis of strangles and the isolation of atypical *Streptococcus equi*. *Vet Rec* 133:215–216.
- Gronbaek L, Angen O, Vigre H, Olsen S: 2006, Evaluation of a nested PCR test and bacterial culture of swabs from the nasal passages and from abscesses in relation to diagnosis of *Streptococcus equi* infection (strangles). *Equine Vet J* 38:59–63.
- Jeong Y-J, Park K, Kim D-E: 2009, Isothermal DNA amplification in vitro: the helicase-dependent amplification system. *Cell Mol Life Sci* 66:3325–3336.
- Kuwamoto Y, Anzai T, Wada R: 2001, Microplate sugar-fermentation assay distinguishes *Streptococcus equi* from other streptococci of Lancefield's Group C. *J Equine Sci* 12:47–49.
- Newton J, Verheyen K, Talbot N, et al.: 2000, Control of strangles outbreaks by isolation of guttural pouch carriers identified using PCR and culture of *Streptococcus equi*. *Equine Vet J* 32:515–526.
- Newton JR, Wood JLN, Dunn KA, et al.: 1997, Naturally occurring persistent and asymptomatic infection of the guttural pouches of horses with *Streptococcus equi*. *Vet Rec* 140:84–90.
- Newton R, Waller A, King A: 2005, Investigation of suspected adverse reactions following strangles vaccination in horses. *Vet Rec* 156:291–292.
- Pooria G, Mohsen A, Amir G, et al.: 2007, Detection of *Helicobacter pylori* by enzyme-linked immunosorbent assay of thermophilic helicase-dependent isothermal DNA amplification. *Diagn Microbiol Infect Dis* 59:243–249.
- Prescott JF, Timoney JF: 2007, Could we eradicate strangles in equids? *J Am Vet Med Assoc* 231:377–378.
- Rozen S, Skaletsky HJ: 2000, Primer3 on the WWW for general users and for biologist programmers. In: *Bioinformatics methods and protocols: methods in molecular biology*, ed. Krawetz S, Misener S, pp. 365–386. Humana Press, Totowa, NJ.
- Sweeney C, Timoney J, Newton J, Hines M: 2005, *Streptococcus equi* infections in horses: guidelines for treatment, control, and prevention of strangles. *J Vet Intern Med* 19:123–134.
- Tang W, Chow WH, Li Y, et al.: 2010, Nucleic acid assay system for tier II laboratories and moderately complex clinics to detect HIV in low resource settings. *J Infect Dis* 201(Suppl 1): S46–S51.
- Timoney JF: 2004, The pathogenic equine streptococci. *Vet Res* 35:397–409.
- Timoney JF, Artiushin SC: 1997, Detection of *Streptococcus equi* in equine nasal swabs and washes by DNA amplification. *Vet Rec* 141:446–447.
- Timoney JF, Artiushin SC, Boschwitz JS: 1997, Comparison of the sequences and functions of *Streptococcus equi* M-like proteins SeM and SzPSe. *Infect Immun* 65:3600–3605.
- Tiwari R, Qin A, Artiushin S, Timoney JF: 2007, Se18.9, an anti-phagocytic factor H binding protein of *Streptococcus equi*. *Vet Microbiol* 121:105–115.
- Tong Y, Tang W, Kim H-J, et al.: 2008, Development of isothermal TaqMan assays for detection of biothreat organisms. *Biotechniques* 45:543–555.
- Vincent M, Xu Y, Kong H: 2004, Helicase-dependent isothermal DNA amplification. *EMBO Rep* 5:795–800.
- Waller AS, Jolley KA: 2007, Getting a grip on strangles: recent progress towards improved diagnostics and vaccines. *Vet J* 173:492–501.
- Webb K, Jolley KA, Mitchell Z, et al.: 2008, Development of an unambiguous and discriminatory multilocus sequence typing scheme for the *Streptococcus zooepidemicus* group. *Microbiology* 154:3016–3024.